

## Adrenal Chromaffin Cell Cultures

For rat-derived cultures, adrenal glands from 7- to 12-day-old Sprague Dawley rats are dissected in ice-cold Hanks Balanced Salt Solution (HBSS).

1. Animals are decapitated (anaesthetize the animals with Ketamine, Ketaset Fort Dodge NDC-0856-2013-01 if decapitation is not an approved protocol).
2. Decapitate and pin the body belly-down. Spray with 70% ethanol.
3. Cut the skin along the spinal column (it's easier starting from the neck) and pull it out to both sides, using a scissor to separate the skin from underlying tissue. The back of the body is now open.
4. Grab the vertical column with forceps and pull it up while cutting along both sides through the ribs. It's better to cut both sides simultaneously (alternating cuts on each side as you work back). Eventually you will see the diaphragm. At this point, open the scissors and place onto the diaphragm approximately 1/3 of the way up from the bottom. Pull the spine up while holding the diaphragm down. The abdominal cavity should tear open exposing, among other organs, the two kidneys with the adrenal glands on top. These should be readily visible.
5. Remove the glands with fine forceps (curved are the best, pinch off the tissue under the glands and pull up) and put them into ice-cold HBSS (Ca<sup>2+</sup>/Mg<sup>2+</sup> free).
6. Adrenal glands are encased by adipose tissue and a capsule. Remove the capsule (using two fine forceps, try to pull the capsule open like a bag of potato chips, then holding the capsule with one tweezer, use the other to "roll" off the gland). Cut the adrenal glands in half (or thirds depending on size of glands).
7. After several washes with HBSS (use a sterile plastic transfer pipette), the tissue is digested with **collagenase IA** solution in 10ml HBSS, Ca<sup>2+</sup>/Mg<sup>2+</sup>-free, (250-350U/ml, Worthington) for about 30-45 min at 37°C with stirring. Stop the digestion as soon as the solution starts to turn cloudy.
8. The digested tissue is rinsed with HBSS and triturated gently in a solution containing 1% heat-inactivated fetal bovine serum and 0.02% DNase I. For trituration use large bore tech-tips, and if needed, medium bore tech-tips (please refer to our "Ventral Midbrain Cultures" protocol for the instruction on how to prepare the tech-tips).
9. Dissociated cells are centrifuged at 1000 g for 3 min to form a pellet and then resuspended in a culture medium comprised of DMEM, 10% fetal bovine serum, 50 U/ml penicillin, 50 µg/ml of streptomycin, and 2mM Glutamine.
10. The cell suspension is plated onto poly-D-lysine- and laminin-coated glass wells in 50 mm dishes (cells from 5 rat 10day-old pups onto 16 dishes) and, after 2 hr, the dishes are flooded

with the \_culture medium (3ml per dish). Please refer to our “Ventral Midbrain Cultures” protocol for the instruction on how to prepare and coat the dishes .

11. Cells are maintained in a 5% CO<sub>2</sub> incubator at\_ 37°C.

All measurements are conducted between 1 and 7 post-plating\_ days.

**Suggestions for plating density for rat and mouse CC cultures:**

3 10-day-old rat pups – 8 dishes  
5 10-day-old rat pups – 16 dishes

2 adult mice - 12 dishes  
10 10-day-old mouse pups - 12 dishes

**CC Media 200ml:**

2ml L-Glutamine  
240ul Pen-Strep  
20ml Fetal Bovine Serum, heat-inactivated  
180ml DMEM

**For 10ml trituration solution:**

10ml HBSS  
30ul DNase stock (final concentration 0.02%)  
100ul Fetal Bovine Serum, heat-inactivated

**Catalog #'s:**

L-Glutamine – 200mM, Sigma # G2150  
Pen-Strep - Sigma # P0781  
Fetal Bovine Serum, heat-inactivated - Gibco # 16140-063  
DMEM - Sigma # D5546  
HBSS – Gibco # 14175-095  
Collagenase Type I – Worthington # CLS-1  
DNase I – Worthington # 2006

**Preparation of DNase I stock solution:**

Reconstitute with HBSS to a concentration of 2,000 U/ml  
(for example: a vial with 20 mg and 3,364 U/mg is reconstituted with 33.64 ml HBSS).  
Store as 500ul aliquots at -80°C.